

Data Acquisition

- Obtain volumetric imaging of the cochlear anatomy
- Determine material properties (conductivity and permittivity)
- Estimate electric loads (shape and amplitude of input current)



Image Processing and Segmentation

- Filter and/or crop scan data
- Identify and segment tissues by type



Solid Modelling

- Convert segmented image stack into an *in silico* replica
- Combine with CAD model of the intracochlear electrode array



Numerical Computation

- Generate appropriate volume mesh
- Assign material properties to corresponding tissues
- Apply electric loads and impose boundary conditions
- Solve for relevant quantities



Postprocessing

- Export and consolidate numerical results
- Calculate derived quantities
- Visualise data using charts, graphics, animations, etc.
- Interpret data in the context of cochlear implant simulations, implant design, and patient outcomes



Validation

- Test volume mesh for convergence
- Perform sensitivity studies on input parameters
- Obtain independent *in vivo* measurements of intracochlear voltages
- Compare *in silico* predictions with *in vivo* measurements